



# Rodenstock Perfalt Tint and Colour Overview.


## Filter Tints




Brown 10 %




Blau 10%




Grey 10%



Brunal  
(brown) 12 %




Rogal  
(rosé) 12 %




Ferral  
(blue grey) 12 %


## Standard Uni Colours




Brown 15 %




Brown 25 %




Brown 50 %




Brown 65 %




Brown 75 %




Brown 85 %




Grey 15 %




Grey 20 %




Grey 50 %




Grey 65 %



Grey 75 %




Grey 85 %




Green 75 %




Green 85 %




Blue 15 %




Blue 25 %




Blue 35 %




Blue 65 %



Yellow 15 %




Amethyst 20 %




Sand/Orange 15 %


## Special Uni Colours




Brown  
Nature 75 %



Grey  
City 75 %




Green  
Pilot 85 %




Fashion  
Black 90 %


## SunContrast




SunContrast  
Amber 65 %




SunContrast  
Amber 80 %




SunContrast  
Brown 75 %




SunContrast  
Brown 85 %




SunContrast  
Green 85 %



SunContrast  
Orange 40 %




SunContrast  
Drive 75 %




SunContrast  
Drive Graduell  
10-75 %

## Mirror Coatings




Silver Moon  
2 %  
combinable with all colours


## Graduell




Brown 10-25%




Brown 10-50%




Brown 10-75%




Grey 10-25%




Grey 10-50%



Grey 10-75%




Blue 10-25%




Blue 10-35%


## Pastella




Sand




Coral




River




Maron



Turquoise



Lagune



Ocean


## Skyline




Black Blue




Black Camel



Terra




Azur



Smoke


## Special Tints



L400 12%



L660 80%



L660 90%

## Coatings for Perfalit Tints

	Coating for Perfalit 1.5			Coating for Perfalit 1.6		
	Silver Moon <sup>1)</sup>	Duralux*	Solitaire	Silver Moon <sup>2)</sup>	Duralux*	Solitaire TopCoat <sup>3)</sup>
<b>Filter Tints</b>	●		●	●		●
<b>Uni Colours</b>	●	●	●	●		●
<b>Special Uni Colours</b>	●		●	●		●
<b>SunContrast</b>	●		●	●		●
<b>Graduell</b>	●	●	●	●		●
<b>Pastella</b>	●		●	●		●
<b>Skyline</b>	●		●	●		●
<b>Special Tints</b>						
- L400			●			●
- L660*			●			

\* Not available for Impression.

We do not recommend to mix two different materials (e.g. 1.5 and 1.6) due to AR matching problems.

<sup>1)</sup> Silver Moon incl. Back Side AR

<sup>2)</sup> Silver Moon incl. Duralux and Back Side AR Solitaire

<sup>3)</sup> also available with Solitaire

Colours may vary slightly.