


# Multigressiv MyView®

Setting a new standard in individualized progressive lenses



**Standard Values for Ordering** - PD, Seg HT, and only now NEAR PD

**Quick Turnaround Time** - 7 working days from the date lenses are ordered from Rodenstock and returned to your lab. GUARANTEED!

**Three progressive zone lengths** - More zone lengths than any other progressive  
 14mm for 14-16mm fitting heights  
 16mm for 16-18mm fitting heights  
 18mm for 18mm and above

MATERIAL	POWER	CYLINDERS	ADD RANGE
MyView® 1.50	-8.00 to +8.00	up to -4.00	0.75 to 3.50
MyView® 1.54 Colormatic	-8.00 to +8.00	up to -4.00	0.75 to 3.50
MyView® 1.60	-10.00 to +8.00	up to -4.00	0.75 to 3.50
MyView® 1.60 Colormatic	-10.00 to +8.00	up to -4.00	0.75 to 3.50
MyView® 1.67	-10.00 to +8.00	up to -6.00	0.75 to 3.50
MyView® 1.67 Colormatic	-10.00 to +8.00	up to -6.00	0.75 to 3.50
MyView® 1.74	-10.00 to +8.00	up to -4.00	0.75 to 3.50

**Rodenstock - German Eyewear.  
 First choice in lenses and frames.**

